

UPPSALA UNIVERSITET

Isomeric fission yield measurements at IGISOL

Outline

- Motivation
- Experimental Facility
- Preliminary Results
- Future Plans



UPPSALA UNIVERSITET

V. Rakopoulos¹ M. Lantz¹ A. Solders¹ A. Mattera¹ A. Al-Adili¹ S. Pomp¹





- A. Jokinen²
- H. Penttilä²
- S. Rinta-Antila²
- V. Simutkin²



JYVÄSKYLÄN YLIOPISTO UNIVERSITY OF JYVÄSKYLÄ



UPPSALA UNIVERSITET FISSION YIELDS

High quality data on independent fission yields

- (I) In Fundamental Physics for:
 - a better understanding of the fission process (model development)
 - studies of neutron rich nuclei far from the line of stability

(II) In Nuclear Energy Applications:

- decrease the uncertainty about spent fuel composition (decay heat, radio-toxicity, storage, reprocessing, transmutation)
- relevant for safety measures

(fission gas production, delayed neutrons, criticality, etc.)

improve burn-up predictions



UPPSALA UNIVERSITET ISOMERIC Yields

Why to bother with the isomers?

- Used to deduce the fission fragment angular momentum.
- Simulations of the astrophysical r-process.
- Simulations of the neutronics and decay heat of nuclear reactors.

e.g: The beta delayed neutron probability from the isomeric state can be an order of magnitude different from that of the ground state (e.g. 0.33% for ⁹⁸Y, 3.5% for ⁹⁸mY).



UPPSALA Motivation

Accurate Fission data FOr Nuclear Safety (AIFONS) project aims at high precision measurements of fission yields

<u>Idea</u>:

To combine fission ion guide with mass identification of ions using a Penning Trap.

Fission Ion Guide is independent to chemistry => ions of every element can be produced.

To identify every nuclide (Z,A) and not only mass selection which can be achieved by the mass separators.









UNIVERSITET IGISOL -JYFLTRAP*





Fission Chamber UPPSALA UNIVERSITET and Ion Guide

- Helium gas jet
 - Stops fission fragments (< 2 MeV)*
 - Recombines fragments to singly charged
 - Guides fragments to SPIG
- SPIG
 - Sextupole ion guide
 - Centres fragments on axis
 - Accelerates by static potentials
- Extractor

Accelerates fragments in steps to 30 kV

*A. Al-Adili et al: "Ion counting efficiencies at the IGISOL facility" arXiv:1409.0714v1







UPPSALA UNIVERSITET IGISOL -JYFLTRAP

























UPPSALA Mass Spectra





UNIVERSITET Proton-induced fission on ²³⁸U

Which nuclides can be separated with the purification trap:

Nuclide	⁸¹ Ge	96 Y	97 Y	⁹⁷ Nb	¹²⁸ Sn	¹³⁰ Sn
T _{1/2}	7.6 s	5.34 s	3.75 s	72.1 m	59.07 m	3.72 m
Jp	9/2+	0-	1/2-	9/2+	0+	0+
lsomer	^{81m} Ge	^{96m} Y	97m Y	^{97m} Nb	^{128m} Sn	^{130m} Sn
Jp	1/2+	8+	9/2+	1/2-	7-	7-
T _{1/2}	7.6 s	9.6 s	1.17 s	58.7 s	6.5 s	1.7 m
E _{level} (keV)	679	1140	667	743	2091	1946



UNIVERSITET Preliminary Results*

Isomeric pair experiment	⁸¹ Ge	96Y	97 Y	⁹⁷ Nb	¹²⁸ Sn	¹³⁰ Sn
2013-06-04	0.06±0.05	1.09±0.11	Unclear		Unclear	0.39±0.27
2013-08-29	0.09±0.05	1.25±0.22	1.85±0.23	0.20±0.60	Unclear	
2014-05-27	0.58±0.27				0.90±0.12	0.66±0.38

*courtesy of D. Gorelov



UPPSALA Future Plans

- Neutron induced independent fission yields of:
- Various actinides (²³⁸U, ²³⁵U, ²³²Th, ...)
- Energy dependence
- Isomeric yield ratios

with high neutron energies (>14 MeV)

Characterisation of the neutron flux on site (winter-spring 2015)







UPPSALA UNIVERSITET





UPPSALA UNIVERSITET

Thank you for your attention!

Financiers for the Uppsala AIFONS project



This work was supported by the European Commission within the Seventh Framework Program through Fission-2010-ERINDA (project no.269499). Strål säkerhets myndigheten swedish Radiation Safety Authority

The Swedish Radiation Safety Authority (SSM)



Svensk Kärnbränslehantering AB Swedish Nuclear Fuel and Waste Management

Co. (SKB)

Swedish Nuclear Physics Meeting 2014, Uppsala



UNIVERSITET Timing of the experiment



TIME

The selected purification cycle length t_p dictates the timing of the experiment. While an ion bunch is purified in the Penning trap, the next bunch is prepared in the RFQ.